Deutsch

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English

INSTRUCTIONS FOR USE

This instrument has been designed for use in combination with the MyLunos Duo ultrasonic scaler and in accordance with the intended use.

Using the specific instruments, the following treatments can be performed:

- Scaling: all procedures for the removal of bacterial plaque deposits and supragingival, subgingival and interdental calculus and for the removal of stains;

- Periodontology: in periodontal treatment for scaling and root-planing/debridement, including cleaning and irrigation of the periodontal pocket.

The instrument must be used only by specialised and properly trained personnel, such as a dentist and/or dental hygienist, able-bodied adults of any weight, age, height, gender and nationality. This instrument must be used in accordance with the correct settings reported in Table 1. The torque wrench supplied with the instrument must only be used to secure and remove the instruments from the MyLunos Duo ultrasonic scaler.

With regard to the intended patient group, patient selection criteria, indications for use, contraindications and warnings, refer to the use and maintenance manual supplied with the ultrasonic piezo-electric scaler device.

The clinical benefits are reported on Chapter 1.2 «Clinical Benefits» of the device's Use and maintenance manual (REF 20441100016).

The terms **WARNING**, **CAUTION** and **NOTE** identify conditions that require particular care and must, therefore, be carefully examined. Specifically:

⚠ **WARNING:** Identifies conditions and practices that present a risk of death or serious injury for the patient and/or user.

⚠ **CAUTION:** Identifies conditions and practices that may cause minor injuries and/or damage to the devices.

NOTE: A condition not associated with a risk or imminent danger.

ⓘ **CAUTION:** The information contained in this document is limited to that strictly necessary for the first use of the product and must be consulted as a supplement to the Use and maintenance manual supplied with the device. Read the Use and maintenance manual supplied with the device carefully, paying particular attention to the «Safety Provisions» chapter, before performing any type of operation on the system and commencing any type of treatment.

⚠ **WARNING:** Before first use and subsequent uses.

The instrument and torque wrench are supplied non-sterile, therefore, before use, they must be reconditioned in accordance with the procedures summarily described below and reported in detail in the «Reprocessing» section of the manual supplied with the device. After each use, the product must be cleaned and sterilised.

⚠ **WARNING:** If an adverse event and/or serious accident attributable to the device occurs during correct and intended use, it is recommended to report it to the Competent Authority and to the manufacturer indicated on the product label.

1 - Cleaning Procedures

IMPORTANT: The cleaning procedures must be performed immediately after each use. Submerge the instrument and/or the torque wrench in demineralised water or the cleaning solution (2% (v/v) ID 213 solution, Dürr Dental) immediately after use.

Do not leave residue or blood deposits on the instruments and the torque wrench, eliminate larger impurities with a disposable cloth or paper towel.

The cleaning procedure can be performed using two different methods: manual cleaning and automatic cleaning. Completely wipe down the external surfaces with FD 322 Premium Wipes for at least one minute until they are visually clean. Make sure that the surfaces are sufficiently moistened.

Note the action time of the cleaning agent.

These two methods are alternative and only one or the other must be performed.

Repeated reconditioning has a minimal effect on these devices and instruments. The end of service life is generally determined by wear or damage due to use.

ⓘ **CAUTION:** Do not use hydrogen peroxide to clean the instruments.

⚠ **WARNING:** To avoid injury to operators, treat the instruments separately.

Particular attention must be paid while handling sharp and potentially infected instruments to avoid the risk of wounds and infections. Health personnel in charge of use and cleaning of the devices must use suitable protection and always use care when handling instruments with sharp edges and pointed tips.

1.1 - Manual Cleaning

NOTE: The instrument and torque wrench must be processed separately. Do not leave the instrument inserted in the torque wrench, do not lay the instrument and torque wrench in the same container, do not use the same brush for the instrument and the torque wrench.

NOTE: Process validated by independent body with a cleaning solution (2% (v/v) ID 213 solution). Use a disposable syringe to aspirate and inject 20 ml of cleaning solution (2% (v/v) ID 213 solution) into Lunos® scaler

instrument's lumen (internal channel) and instruments torque wrench hard-to-reach areas (through holes, internal cavities, grooves and crevices).

Remove the syringe.

Afterwards place the individual parts in a cleaning and disinfecting bath (non-fixing, aldehyde-free) for the stated reaction time (for 2% ID 213, min. 5 minutes, max. 2 hours). Make sure that all parts are covered.

In the disinfectant bath, brush all exterior and interior surfaces completely with a soft hygienic brush.

Comply with the action times of the cleaning agent and disinfectant (minimum 1 minute).

NOTE: Thoroughly brush, for minimum 1 minute each, all the following parts of the instruments and torque wrench:

- Through-holes and internal channels;
- Hard-to-clean areas such as sharp edges and in particular, the interstices between the cutting cusps;
- External metal ring;
- Internal cavities, grooves and fissures.

After the action time prescribed by the manufacturer of the cleaning and disinfectant solution, rinse all components under water for minimum 1 minute (temperature < 35 °C).

Use a disposable 20 ml syringe to aspirate and inject tap water for minimum 1 minute (temperature < 35 °C) into Lunos® scaler instrument's lumen (internal channel) and instruments torque wrench hard-to-reach areas (through holes, internal cavities, grooves and crevices).

If the surfaces are not visibly clean, repeat cleaning and inspect again. Eliminate damaged instruments.

To complete the manual cleaning use a disposable syringe to aspirate and inject 20 ml of disinfectant solution (2% (v/v) ID 213 solution) into Lunos® scaler instrument's lumen (internal channel) and into instruments torque wrench hard-to-reach areas (through holes, internal cavities, grooves and crevices).

Remove the syringe.

Afterwards place the individual parts in a cleaning and disinfecting bath (non-fixing, aldehyde-free) for the stated reaction time (for 2% ID 213, min. 5 minutes, max. 2 hours). Make sure that all parts are covered.

After the action time prescribed by the manufacturer of the cleaning and disinfectant solution, rinse all components under deionised water for minimum 1 minute (temperature < 35 °C).

Use a disposable 20 ml syringe to aspirate and inject deionised water for minimum 1 minute into Lunos® scaler instrument's lumen (internal channel) and into instruments torque wrench hard-to-reach areas (through holes, internal cavities, grooves and crevices).

1.2 - Automatic Cleaning

NOTE: The instrument and torque wrench must be processed separately. Do not leave the instrument inserted in the torque wrench.

Automatic cleaning involves the use of a thermal disinfectant and the following materials:

- Alkaline detergent: neodisher® FA (0.2 % v/v);
- Neutralising liquid: neodisher® Z (0.1 % v/v);
- Metal basket;
- Adaptors.

NOTE: Make sure that the accessories are appropriately blocked in the basket and cannot move during washing. Any blows could damage them. Position the instruments in a way that the water can flow through all the surfaces, even internal.

⚠ WARNING: Avoid overloading the thermal disinfectant as this could compromise cleaning effectiveness.

1 Place the accessories in a metal basket or use the special adaptors (supplied as optionals).

2 Set the sequence and parameters for the cycle in the thermal disinfectant:

- 1 min, Rinse with cold water;
- 5 min, Wash with alkaline detergent at 55 °C ±2 °C;
- 1 min, Neutralisation with neutralising liquid and deionized water at 32 °C ±2 °C;
- 1 min, Rinse with deionized water at 32 °C ±2 °C;
- 5 min, Thermal disinfection at 93 °C with deionized water.

Thermal disinfection has not been tested experimentally. In compliance with ISO 15883-1, Table B.1 [4] thermal disinfection at a temperature of 90°C for 5 min determines a value of A0 3000.

⚠ CAUTION: The use of alkaline detergents removes the lubricant from the torque wrench, causing malfunctions and increasing its wear. Therefore, it is important to lubricate the torque wrench after each cleaning procedure using medical-grade lubricants (see Paragraph 1.3). The dilution of the alkaline detergent required for automatic cleaning is such that it does not damage the silicone parts of the torque wrench (o-ring).

1.3 - Lubrication

Before sterilisation, the torque wrench must be lubricated with a commercial medical-grade lubricant.

NOTE: In U.S. , adopt a medical-grade lubricant approved by FDA, for example, the NSK Pana-Spray or similar.

⚠ CAUTION: Do not lubricate the instruments.

- The lubricant must be applied by spraying it directly onto the peripheral contact surface inside the torque wrench, as indicated in Figure 2;
- After having applied the lubricant, remove any excess oil using a clean lint-free cloth.

1.4 - Packaging

The instrument and torque wrench must be sterilised using a standard, appropriately sized medical steam sterilisation pouch.

⚠ CAUTION: The instruments and torque wrench must be packaged separately. Do not package multiple instruments and/or torque wrenches in the same pouch.

⚠ CAUTION: Make sure the pouch is big enough to hold the instrument without tensioning the seals and without ripping the packaging.

1.5 - Sterilization

Sterilise only by pre-vacuum steam autoclaving.

The steam sterilisation parameters have been validated by the manufacturer to guarantee a sterility level (Sterility Assurance Level - SAL) of 10⁻⁶.

⚠ CAUTION: If multiple instruments need to be sterilised in the autoclave in a single cycle, do not exceed the maximum allowed load.

⚠ CAUTION: Risk of contamination. Do not use gravity autoclaves to sterilise the instruments. The operating cycle of gravity autoclaves does not guarantee the suitable sterilisation of the internal channel, of the cavities and hard-to-reach points.

⚠ CAUTION: Use only a pre-vacuum autoclave to sterilise the instruments and the torque wrench. Do not use other methods of sterilisation insofar as they may be incompatible with the materials used to produce the instruments.

Do not use the following methods of sterilisation: ethylene oxide sterilisation, hot air sterilisation, flash autoclaving, STERRAD sterilisation, sterilisation with STERIS or similar sterilisation systems. Do not use the following substances/systems to sterilise the instruments: hydrogen peroxide, peracetic acid, Formaldehyde, Glutaraldehyde based systems or other equivalent solutions/systems.

1.5.1 - Sterilisation - Parameters

Minimum sterilisation parameters to guarantee a sterility level (Sterility Assurance Level - SAL) of 10⁻⁶:

⚠ WARNING: these sterilisation parameters, validated by an independent laboratory, refer exclusively to a PRE-VACUUM autoclave operating cycle.

- Type of cycle: PRE-VACUUM - (3 pre-vacuum stages - 60mbar).
- Minimum temperature: 132 °C + 3 °C (270 °F + 5 °F).
- Pressure: 2.9 bar (42 psi).
- Minimum exposure time (in package): 4 minutes.
- Minimum drying time: 20 minutes.

⚠ CAUTION: Do not exceed a temperature of 135 °C.

All of the stages of sterilization must be performed by the operator in compliance with ISO 17665-1:2007, ISO 556-1:2002 and ANSI/AAMI ST 46:2002.

2 - Important information relative to safety of instruments

⚠ WARNING: Breakage and wear of instrument.

On rare occasions, the oscillations at high frequency and wear may cause an instrument to break. Do not bend, reshape or sharpen the instruments in any way. The alteration of an instrument may cause it to break. Altered instruments must never be used. If an instrument breaks during use, recover and remove all fragments of the instrument from the patient, with maximum care. During the treatment, ask the patient to breathe through their nose and/or use a dental dam to avoid the risk of swallowing any fragments which may be generated in case of breakage of the instrument.

⚠ WARNING:

• The instrument is a consumable. Thoroughly inspect the instrument before and during every treatment, looking for damage and/or excessive wear. Do not use an instrument if it is scratched or rusted. The instrument may break during use. If damage or deterioration in performance is noted, replace the instrument with a new one. Use of a worn instrument shall lower performance in terms of cutting power.

NOTE: A deformed instrument may preclude the complete functionality of the device used, and, depending on the model, a message may be displayed.

• Before the treatment, check that the instrument is firmly secured to the handpiece. Properly secure the instrument onto the handpiece using the torque wrench.

• Do not modify the shape of the instrument in any way. Any bending or prying of the instrument may cause it to fracture. Never use deformed instruments.

• Do not attempt to sharpen used instruments insofar as they may break during use.

• Replace the instrument/s only with original spare parts. Use of instruments other than original instruments shall void the warranty of the device. Never use instruments other than original instruments insofar as they will damage the device and may cause injury to operators or the patient. The use of instruments other than original instruments shall damage the threaded pin of the handpiece, thus compromising the secure attachment of instruments, even if original, to the handpiece.

• It is advisable to avoid the application of excessive force or prolonged contact of the instrument on the soft tissues to avoid thermal damage and / or injury.

• Let the ultrasonic vibrations work, do not exert excessive pressure on the instruments during use. Apply a light force on the instrument to obtain the best efficiency.

• The instrument must always be kept moving. If the instrument is blocked, it can cause overheating of the treated part. It is recommended to use continuous movement to minimize contact between the tip and the treated surface. Do not block it against the tissue so as not to cause it to overheat. It is advisable to use high levels of irrigation as the power level increases.

• The instruments vibrate with a longitudinal oscillation, with forward and backward movement. During treatment, always keep the instrument tangential to the tooth surface. Move the handpiece back and forth while applying light lateral pressure.

• Do not aim the instrument directly on the surface of the enamel or implant. Position the tip/operative part only tangentially to the surface of the tooth or implant.

• When the instrument is used in the interproximal spaces, do not block the instrument or lever with the working tip. The instruments must be left free to vibrate.

• Check the threaded parts of the instrument and those on the handpiece. These parts must be thoroughly cleaned.

⚠ WARNING: Hospital waste

Worn or damaged instruments must be treated as hospital waste.

3 - Use of the torque wrench.

⚠ WARNING: During operations to secure and remove the instrument, the user must pay special attention to the sharp and cutting edges of the instruments.

⚠ WARNING: Hospital waste.

A worn or damaged torque wrench must be treated as hospital waste.

3.1 - Securing the instrument

For correct use of the torque wrench, operate as follows:

- 1 Insert the pre-selected instrument all the way onto the handpiece;
- 2 Insert the instrument into the torque wrench, as illustrated in Figure 1;
- 3 Firmly grip the body of the handpiece.

⚠ CAUTION: Do not grip the handpiece by the end and/or cord, but only by the body (Figure 1 - Ref. C). Do not rotate the body while securing the instrument. Firmly grip the body of the handpiece and turn the torque wrench only.

4 Turn the torque wrench in a clockwise direction until the friction clicks (the external body of the torque wrench turns with respect to the body of the handpiece, emitting mechanical "CLICK" sounds).

5 The instrument is now properly locked in.

3.2 - Removing the instrument

Position the torque wrench on the instrument and screw in an anti-clockwise direction.

MyLunos Duo	ULTRASOUND Funktionen / functions	
	Scaler	Perio Scaler
SCALING-INSTRUMENTE / SCALING INSTRUMENTS		
S1-S	1 - 3	
S2	1 - 5	
PERIO-INSTRUMENTE / PERIO INSTRUMENTS		
P10 - P11 - P12	1 - 4	1 - 3
P16R - P16L		1 - 3

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Mectron S.p.A.